



# International Testing Methods for food relevant microorganisms

Wolfgang E. Schmidt , MERCK KGaA , Darmstadt , Germany

# Standard Methods



**ISO**




**FDA / BAM**

- ❖ Salmonella
- ❖ Listeria
- ❖ Campylobacter
- ❖ pathogenic E.coli
- ❖ Cronobacter

Bacillus cereus  
Clostridium perfringens  
Coliforms / E.coli  
Enterococci  
Legionella

# Salmonella

- Rod shaped
- Gram negative enterobacteria
- Salmonella species are motile and produces hydrogen sulfide
- Usually not able to ferment Lactose
- Near related to Escherichia coli
- Can be found world wide



**DANIEL ELMER SALMON,**  
**D.V.M. (1850-1914)**

# Salmonella

- Leading cause of food borne bacterial diseases in many countries
- ~2500 Serovars
- Divided in 2 species:
  - *S. enterica* - non-typhoid serovars
    - 6 Subspecies
  - *S. bongori* - typhoid serovars



# Salmonella



- High environmental resistance (can survive 73 days at 2°C or at pH 3.5)
- Deep freezing don't kill Salmonella
- Infections result from contaminated meat (poultry!!!), eggs, milk products (cheese, cream, ice cream, etc.), spices, etc.

## Major outbreaks of Salmonellosis:

Year	Country	Food	Serovar	No. Infection	No. Deaths
1984	Canada	Cheddar Cheese	<i>S. typhimurium</i>	2,700	0
1985	USA	Pasteur. Milk	<i>S. typhimurium</i>	16,284	11
1988	Japan	Cooked eggs	<i>Salmonella spp.</i>	10,476	N.S.
1993	Germany	Paprika Chips	<i>S. saint-paul</i>	1,000	0
2007	Germany	Cream	./.	over 260	8

# Salmonella

## Enrichment culture

- Buffered peptone water
- Tetrathionate broth
- Rappaport Vassiliadis



# DIRECT OR PRE-ENRICHMENT ?




TYPE OF SAMPLE	RECOVERY % ENRICHMENT		REFERENCE
	PRE-	DIRECT	
<b>RAW MEAT</b> (1983)	<b>83</b>	<b>56</b>	<b>EDEL &amp; KAMPELMACHER</b>
	<b>90</b>	<b>58</b>	<b>VASSILIADIS ET AL. (1981)</b>
	<b>91</b>	<b>62</b>	<b>GABIS &amp; SILLIKER (1974)</b>
	<b>60</b>	<b>55</b>	<b>EDEL ET AL. (1976)</b>
<b>CHICKEN</b>	<b>96</b>	<b>10</b>	<b>VASSILIADIS ET AL. (1972)</b>
<b>CARCASSES</b>	<b>43</b>	<b>47</b>	<b>COX ET AL (1978)</b>
	<b>99</b>	<b>99</b>	<b>COX ET AL. (1979)</b>
<b>ANIMAL FEED</b>	<b>43</b>	<b>36</b>	<b>JUVEN ET AL. (1984)</b>
	<b>40</b>	<b>33</b>	



# Buffering Capacity of Buffered Peptone Water

Laboratorio de Sanidad Animal, Barcelona.  
Departamento de Microbiología. M.Saco.



**MERCK**


# Salmonella

## Isolation media


- XLD Agar
- Rambach
- Any other agar



# ISO Standard 6579 for Detection of Salmonella




**Discover the new Readybag pouches  
with dehydrated culture media for  
fast , flexible and easy pathogen testing**

An abstract graphic consisting of several overlapping, curved, light blue and white shapes that sweep upwards and to the right, creating a sense of motion and modernity.






# Readybag concept

## Speed up and simplify your workflow!

- Pre-weighed
- Gamma-irradiated
- No media preparation needed!
- Eliminates autoclavation
- Sterile lab water provided with ELIX advantage



# Pathogen enrichment - Traditional workflow

1	 	<h2>1. Media preparation</h2> <p>Weighing + Dissolving media</p> <p>5 min</p>
2	 	<ul style="list-style-type: none"> <li>•Add water</li> <li>•Stir and dissolve</li> <li>•Load autoclave</li> <li>•Autoclaving</li> <li>•Cool down autoclave</li> <li>•Unload autoclave</li> <li>•Cool down media</li> </ul> <p>01 min 05 min 02 min 25 min 30 min 05 min 30 min</p> <p>98 min</p>
3		<p>Supplementation Selective Media</p> <p>6 min</p>

# Pathogen enrichment - Traditional workflow


4



## 2. Food sample handling


Weighing food sample 3 min

5



Adding media 3 min

6






Homogenization

2 min

-----  
**117 min**

# Pathogen enrichment – Readybag Workflow

1		Weighing food sample	3 min
2		Add Readybag content	1 min
		Add sterile water	2 min
3		Homogenization	2 min ----- <b>8 min</b>



# High throughput sample solution – ELIX + Heated Water + Readybag Workflow

1



Weighing food sample

2



1. ELIX with tank and Pump
2. Commercial available heater to heat water to approx. 35-37 C
3. Dispenser for eg 225 or 3375 ml water


3



## Homogenization


1. Add Readybag to food sample
2. Add 225 ml or 3375 ml water
3. Blend sample
4. Incubate

# FDA-BAM Standard for Detection of Salmonella



# Immunological Rapid Test Method -


Singlepath® Salmonella for the Detection of Salmonella in Food



Salmosyst is a 2-step enrichment system designed for sublethally damaged *Salmonella*.

Salmosyst is TT-Broth divided in :

- a.) Non-selective Broth
- b.) Selective Enrichment  
as slow-release tablets



## a.) Non-selective Broth

<b>Tryptone</b>	<b>5.0 g</b>
<b>Meat-Peptone</b>	<b>5.0 g</b>
<b>NaCl</b>	<b>5.0 g</b>
<b>Calciumcarbonate</b>	<b>10.0 g</b>

## b.) Selective Enrichment as slow-release tablet

<b>Potassium-tetrathionate</b>	<b>0.2</b>
<b>Oxbile</b>	<b>0.08</b>
<b>Brillantgreen</b>	<b>0.0007</b>
<b>Calciumcarbonate</b>	<b>0.1</b>



# Salmosyst ® 24 h Method



Day  
1



25 g/ml TEST SAMPLE  
IN SALMOSYST® BROTH  
6 h / 35 - 37 °C

ADD 1 SALMOSYST® SELECTIVE  
TABLET TO 10 ml SALMOSYST®  
BROTH AND CONTINUE INCUBATION  
18 h / 35 - 37 °C



1  
Medium  
only



Day  
2

SUBCULTURE ON  
SALMONELLA AGAR'S



# Salmonella Detection from Naturally Contaminated Food Samples by ISO Method in Comparison to Salmosyst® - Rambach® Agar Method



Food Samples	Number of Samples	No. of positive samples (%)	
		ISO Method	Salmosyst® - Rambach® Method
Poultry	85	18 (21.2)	26 (30.5)
Meat	246	14 (5.7)	14 (5.7)
Eggs	188	7 (3.7)	7 (3.7)
Sensitivity:		81.2 %	97.9 %
Total Analysis Time:		4-6 days	2 days
No. of media:		5	2

Study published in Applied and Environmental Microbiology,  
May 1995, p. 1996 - 1999 by Giuseppe Giammanco et al.

## Detection of Salmonella using Rambach® Agar

**Traditional Salmonella media are based on  
 $H_2S$  reaction:**

**Disadvantage:**


**Citrobacter and Proteus mimic Salmonella**

**High number of false - positives**




# **Citrobacter and Proteus mimic Salmonella**

## **High number of false - positives**



**Hektoen**



**XLD**



**BSA**

# SALMONELLA

## RAMBACH AGAR



### Differential system


Red	Salmonella
Violet	Coliforms
Colourless	Other bacteria

### Selective system

Sodium deoxycholate

### Principle


1. Salmonellae dissimilate propylene glycol producing acid. Neutral red gives colonies red colour
2.  $\beta$  galactosidase positive organisms split chromogenic substrate producing blue-violet colonies



# Rambach® Agar



**Clear, Easy and Reliable Identification of Salmonella  
in Food Samples**




**Rambach® Agar # 1.07500**

# Rambach® Agar



## ISO Standard 6579 for Detection of Salmonella




ISO 6579:2002, Amendment 1, Detection of *Salmonella* spp. in animal faeces and in samples from primary production stage

- Animal faeces (such as from poultry pigs cattle)
- Environmental samples in the area of primary production stage (such as dust)

# MSRV Medium



- the detection principle is based on the motility of *Salmonellae* to migrate into the semi-solid medium
- the motility of other organisms is largely inhibited by magnesium chloride, malachite green and Novobiocin and the enhanced incubation temperature of 41.5°C



# Salmonella

## Recommended media


- MSRV-Agar
- DIASALM




# Semi-solid testing Method acc.to ISO 6579



Day  
1



25g/ml TEST SAMPLE  
IN 225 ml BPW  
16 - 20 h / 35 - 37 °C



Day  
2

DIASALM OR MSRV  
24 h / at 41.5 °C



# SALMONELLA DETECTION FROM NATURALLY CONTAMINATED RAW MEAT AND MEAT PRODUCTS BY ISO METHOD IN COMPARISON TO DIASALM METHOD



No. / TYPE of samples	NO. OF POSITIVE SAMPLES	
	MODIFIED ISO (BPW+RVS+XLD)	DIASALM
136 POULTRY	101	134
SENSITIVITY	(74.3 %)	(98.5 %)
102 PORK + BEEF	97	101
SENSITIVITY	(95.1 %)	(99 %)
TOTAL ANALYSIS TIME	3 days	2 days
NO. OF MEDIA	3	2

**STUDY PUBLISHED IN CULTURE MEDIA, DE WARE(N) CHEMICUS  
28 (1998) 34 - 43 BY L. DE ZUTTER *et al.***

# SALMONELLA DETECTION FROM NATURALLY CONTAMINATED RAW MEAT AND MEAT PRODUCTS BY ISO METHOD IN COMPARISON TO MSRV METHOD



No./TYPE of Samples	NO. OF POSITIVE SAMPLES		
	ISO + MSRV	ISO	MSRV
913 (80 % CHICKEN MEAT)	309 (33.8 %)	251 (27.5 %)	291 (31.8 %)
SENSITIVITY		81.2 %	<b>94.5 %</b>
TOTAL ANALYSIS TIME		4 - 6 days	2 days
NO. OF MEDIA		5	2

**STUDY PUBLISHED IN ARCHIV FÜR LEBENSMITTELHYGIENE  
OCTOBER 1993, P. 125 - 127 BY P. PLESS et al.**

# XLT4-AGAR



## Differential system

Black	Salmonellae
Black with yellow halo	Citrobacter
Colourless with yellow halo	Xylose-Lactose-Sucrose positive bacteria
Colourless	Other bacteria e.g. Shigella


## Selective system

Tergitol 4 inhibites Proteus

## Principle

H<sub>2</sub>S formation via thiosulphate  
and iron produces black colonies

Dissimilation of sugars produces acid resulting in phenol red shift of red to yellow



## XLT4 Agar:

### An Improved Selective Agar Medium for the Isolation of Salmonella

R.C. Miller, USDA, 1991 Poultry Science 70



Salmonella isolation rates observed on each medium examined from 91 Salmonella-positive drag-swab samples collected from 7 table egg layer flocks.

Variable	Medium				
	XLT4	XLD	XLDN <sup>1</sup>	BGA	BGAN <sup>1</sup>
Number positive out of 91	89	27	76	65	77
Percentage positive	98 %	30 %	84%	71 %	85 %
Average purity of Salmonella colonies <sup>2</sup>	3+ to 4+	1+	2+	1+ to 2+	2+ to 3+

<sup>1</sup> = Medium supplemented with Novobiocin

<sup>2</sup> = 1+ = 1 to 25 % were Salmonella  
= 2+ = 25 to 50 % were Salmonella  
= 3+ = 51 to 75 % were Salmonella  
= 4+ = 76 to 100 % were Salmonella

# XLD-Agar , XLT 4 Agar and Rambach Agar

Interpretation of colony colours :

	Colony color		
Organisms	XLD Agar	XLT4 Agar	Rambach Agar
<i>Citrobacter freundii</i>			
<i>Proteus</i>			
<i>Salmonella</i>			



# Detection Methods for Listeria and L.monocytogenes



# Established Pathogen: *Listeria*


- Genus *Listeria* comprises 6 different species
- Only *L. monocytogenes* can be pathogenic for humans and animals
- **Disease:**

Systemic Infections:	Meningitis, Encephalitis, Septicemia
Local Infections:	Gastroenteritis
- **Frequency:** Only 7 / 1 Mio. individuals suffer from Listeriosis  
But: Mortality rate is extremely high (up to 25% of patients die)
- **Resistance:** Multiply at 2°C, survive at many preservation methods (10% osmolarity, pH 5.0, 55°C)
- All *Listeria* species serve as indicator for improper hygiene conditions

# Listeria characteristics

Genus Listeria comprises 10 different species

- Listeria monocytogenes      *Listeria marthii*
- Listeria ivanovii      *Listeria rocourtiae*
- Listeria innocua      *Listeria weihenstephanensis*
- Listerai seelegeri      *Listeria fleischmannii*
- Listeria welshimeri
- Listerai grayi





## Listeria disease - Listeriosis

- Listeriosis is an invasive infection caused when contaminated food is ingested.
- Symptoms vary, ranging from flu-like symptoms to meningitis and encephalitis.
- Mortality rate for Listeriosis is 23%.
- The infection usually occurs sporadically.

# Media for enrichment of Listeria



- Half FRASER broth
- FRASER Listeria Enrichment broth
- Listeria Enrichment broth
- Buffered Listeria Enrichment Broth
- UVM Listeria Selective Enrichment broth
- PALCAM Selective Enrichment broth

# Listeria Enrichment Media

- \* Esculin hydrolysis = blackening broth in the presence of Esculin positive organisms like Listeria = is used in UVM , FRASER and PALCAM as presumptive indicator for Listeria spp.




- \* Colour development is depending on level of Listeria present in the sample . To use colour development only is not reliable , therefore all broth cultures must be sub-cultured to plating media . Usually a level of 10.000 c.f.u./ml does not result in blackening of broth

- PALCAM Agar
  - Listeria spp.
- OXFORD Agar
  - Listeria spp.
- Listeria selective agar acc. to AGOSTI & OTTAVIANI
  - Listeria monocytogenes
  - **New mandatory isolation medium in the revised ISO Standard 11290 ( 2004 )**

# Listeria Plating Agars

- \* Common used plating media are PALCAM and OXFORD Agar . Both media are based on the colour development of hydrolysed esculin .




- \* PALCAM Agar contains Mannitol to differentiate Enterococci from Listeria by Mannitol fermentation



# Chromocult® Listeria Selective Agar



**Chromogenic Culture Medium for the Detection and Enumeration of *Listeria monocytogenes* in Foods.**



**Chromocult® Listeria Selektiv Agar # 1.00427**

# Chromocult® Listeria Selective Agar



**Chromocult® Listeria Selective Agar complies with Agar Listeria according to Ottaviani and Agosti (ALOA®) in line with the recommendations of ISO 11290 (2004) and FDA / BAM (2003)**

**1.00427 Chromocult®  
Listeria Selective Agar Base**



**1.00432 Chromocult®  
Listeria Selective Supplement**

**1.00439 Chromocult®  
Listeria Enrichment Supplement**




# ISO 11290-Part 1: Detection Method (*L. monocytogenes* in Food)




Day  
1

Streak out  
on **Chromocult® +  
PALCAM**  
37°C for 24 - 48 h

25 g Test samples  
in 225 ml **1/2 strength FRASER**  
30°C for 24 h



Day  
3 - 4



0.1 ml in 10 ml **FRASER**  
35 or 37°C for 48 h



Day  
3 - 5

**Biochemical Confirmation**

Streak out  
on **Chromocult® +  
PALCAM**  
37°C for 24 - 48 h




**Biochemical Confirmation**

Day  
5 - 7



# ISO 11290-Part 2: Enumeration Method (*L. monocytogenes* in Food)

Day  
1



Day  
2 - 3

Count all colonies showing a  
blue-green colour with opaque halo



Day  
4 - 5

Biochemical Confirmation

# FDA-BAM Standard for the detection of Listeria



Day  
1

25 g/ml TEST SAMPLE  
IN 225 ml **Buffered LEB**  
30 °C for 4 h



Day  
2-3

**After 24 h**

add selective  
agents  
30 °C for 44 h

**After 48 h**


STREAK OUT  
ONTO **Chromocult**®  
37 °C  
FOR 24 - 48h

STREAK OUT  
ONTO **OXFORD**  
35°C  
FOR 24 - 48h

STREAK OUT  
ONTO **Chromocult**®  
37 °C  
FOR 24 - 48h

STREAK OUT  
ONTO **OXFORD**  
35°C  
FOR 24 - 48h

Day  
3-5



PICK 5 CHARACTERISTIC COLONIES AND SUBCULTURE ON TSYE AGAR 35 OR 37°C FOR 24 h

Day  
6

**BIOCHEMICAL CONFIRMATION:** GRAM (+), CATALASE (+), MOTILITY (+), CAMP (+),  
β HAEMOLYSIS (+), GLUCOSE (+), RHAMNOSE (+), XYLOSE (-)

# USDA Standard for testing of Listeria in Red Meat, Poultry, Eggs, Environmentals



Day  
1

25 g/ml TEST SAMPLE  
IN 225 ml **UVM 1**  
30 °C FOR 24 h




Day  
2



0.1 ml to 10 ml to  
**UVM II or FRASER Broth**  
35 °C for 24 h

Day  
3-5

STREAK OUT  
ONTO **OXFORD AGAR**  
35 °C FOR 24-48h



PICK OFF 5 -20 CHARACTERISTIC COLONIES AND SUBCULTURE ON HORSE BLOOD AGAR 35°C FOR 24 h

Day  
6

**BIOCHEMICAL CONFIRMATION:** GRAM (+), CATALASE (+), MOTILITY (+),  
CAMP (+),  $\beta$  HAEMOLYSIS (+), GLUCOSE (+), RHAMNOSE (+) , XYLOSE (-)

# Singlepath® L' mono - Screening

Day  
1

25 g/ml TEST SAMPLE  
IN 225 ml **1/2 STRENGTH FRASER**  
30°C FOR 24h



Day  
2

0,1 ml IN 10 ml **LEB** or  
**Full Fraser** OR **UVM**  
30° / 37°C FOR 21 - 24 h



Day  
3

***L. monocytogenes***  
not present



150µl



***L. monocytogenes***  
present

**Confirmation**



**ALOA®**  
37°C for 24 - 48h

# Listeria Identification Tests

\*

Characteristic	Reaction
Motility (at 20 - 25 °C)	+
Oxygen requirement	facultative
Growth at 35 °C	+
Catalase activity	+
Hydrogen sulphide production	–
Acid from glucose	+
Methyl red reaction	+
Voges-Proskauer reaction	+
Indole production	–
Citrate utilisation	–
Urease activity	–

# Singlepath® L' mono - Confirmation




ALOA® Agar

suspend 1 - 3 presumptive colonies in **250 µl** BHI or CASO or Fraser or bLEB Broth

1 h at 37°C



PALCAM Agar



OXFORD Agar

150 µl



*L. monocytogenes*  
not confirmed




*L. monocytogenes*  
confirmed

**Fastest manual confirmatory test for *Listeria monocytogenes***

# Upcoming Pathogen: *Campylobacter*




- Found in birds, wild and domestic animals
  - Most common bacterial cause of food-borne diseases
  - Thermophilic *C. jejuni* and *C. coli* cause fever and gastroenteritis
  - Duration: 2 - 4 days; incubation period: 1 - 7 days
  - 1/1000 patients develop Guillian Barré Syndrom (Paralysis for weeks - months)
  - Infections result from consumption of contaminated poultry, raw milk products
  - Grows under microaerophilic conditions (difficult to culture in a lab)
- Problem: Many food labs cannot analyse *Campylobacter* in foods because of no knowledge of handling of these pathogens

# Detection of Campylobacter in Food acc. to FDA-BAM and ISO





# Immunological Rapid Test for Detection of Campylobacter in Food by Singlepath® *Campylobacter*

Day 1	only 2 days	<div>25g/ml Test Sample in 225 ml Bolton Broth</div> <div>37°C for 4 h 41,5°C for 44 h</div>		only 1 medium
Day 2				
Day 3 & 4	NO Campylobacter present			YES Campylobacter present ↓
Day 5	If positive, confirmation required by culture and biochemical tests		<div>Streak out on CCDA, 37°C for 48 h</div>	



Enterobacter sakazakii

## ***Enterobacter sakazakii***

- **Considered as a clear species since 1980**
- ***E. sakazakii* is an opportunistic pathogen**
- **Risk for new-borns, especially low-birth weight infants**
- **Causes severe neonatal sepsis, meningitis**
- **Infant formula food has been a vehicle and cause of *E. sakazakii* illness**

# Enterobacter sakazakii


- Powdered infant milk formulas is not a sterile product and may occasional contain pathogen.
- One of the pathogen can be *Enterobacter sakazakii*
- Normally Enterobacter are not harmful and in the past it was not easy to detect this specific Enterobacter



# Chromocult® Enterobacter sakazakii Agar



## ISO TS 22964



# Enterobacter sakazaki Agar



During the past decade there has been a significant increase in the use of chromogenic substrates in isolation media (7).

It has been reported that 100% of *E. Sakazakii* were positive for  $\alpha$ -D-glucosidase and that 100% of other *Enterobacter* species were negative for this enzyme (8).


Based on this observation 5-bromo-4-chloro-3-indolyl- $\alpha$ -D-glucopyranoside (X- $\alpha$ -Glc, 6) and 4-methylumbelliferyl- $\alpha$ -D-glucoside (9) was added to a basal medium to differentiate *E. sakazakii* strains from other members of the *Enterobacteriaceae*.

The enzyme  $\alpha$ -glucosidase hydrolyses X- $\alpha$ -Glc giving blue coloured colonies on DFI agar (Oxoid, UK), ESIA agar (AES, France) or Chromocult ES agar (Merck, Germany), which are commercially available.

# Chromocult® Enterobacter Sakazakii Agar

**Chromogenic Culture Medium for the Detection and Enumeration of *Enterobacter sakazakii* in Milk Powder and Powdered Infant Formula.**

***E. sakazakii***



**Chromocult® Enterobacter Sakazakii Agar # 1.00873**



***Bacillus cereus***



# Bacillus cereus – Sources



**Humans** are not a significant source of food contamination by B.cereus. This organism already exists on many foods and can therefore be transiently carried in the intestine of healthy humans.

**Animals** can carry B.cereus on parts of their body. May cause mastitis in cows.

**Raw Foods** of plant origin are the major source.

The ability of spores to survive dried storage and the thermal resistance of spores, means that most ready-to-eat foods will contain B.cereus.

Strains producing Emetic toxin grow well in rice dishes and other starch foods whereas strains producing diarrhoeal toxin grow in a wide variety of foods from vegetables and salads to meat and casseroles.

Numerous dried herbs and spices and dehydrated foods have been shown to contain B.cereus

# Bacillus cereus – Sources



## Environment

B.cereus is widely distributed in nature and can be found in soil, dust, Air, water and decaying matter. Its ability to form spores allows survival through all stages of food-processing, other than retorting.

## Transmission Routes

Ingestion of contaminated food

## Treatment

Usually no treatment is given. Fluids may be administered when Diarrhoea and vomiting are severe

# Bacillus cereus – Foodborne illness



## 2 Types of Infections:

### 1. **Diarrheal illness** = **Toxicoinfection**

Enterotoxin(s) produced during vegetative growth of *B. cereus* in small intestine

Associated with ingestion of *B. cereus* producing heat-labile toxins  
(occurs within 8-12 hours)

### 2. **Emetic illness (Vomiting)** = **Intoxication**

Cereulide Toxin preformed in food

Usually associated with the ingestion of a heat-stable toxin from contaminated rice  
(occurs within 1-6 hours)



INADEQUATE REFRIGERATION for several HOURS



PREPARATION OF FOODS IN ADVANCE



POOR PERSONAL HYGIENE



INADEQUATE HEATING OR COOKING



PROLONGED USE OF WARMING PLATES  
(PROMOTES Toxine PRODUCTION)

# FOODS AT RISK



FISH

Boiled cod, Lobster pate



MILK

Pasteur.milk, Infant cream



MEAT

Meat loaf, boiled sausage, beef stew, veal broth

CEREALS

Fried rice, cooked rice and pasta dishes

POULTRY

Turkey loaf, roast turkey, barbecued chicken

EGG


Omelette

DESSERTS

Vanilla slices, vanilla pudding, sauce

VEGETABLES

Vegetable sprouts, pea soup, mashed potatoes



HERBS/SPICES

Many foods

# Bacillus cereus: Sources of infections

Meat products, soups, milk & milk products, vegetables, puddings & sauces

– for **diarrhoeal syndrome**



Rice, pasta, noodles, pastry, starchy products

– for **emetic syndrome**



# Outbreaks / infections due to *B. cereus*




## **Most Susceptible Populations**

- Young people (<19 years)
- Elderly (>60 years)

## **Infectious dose:**

- $5 \times 10^4$  to  $10^{11}$  cells / dose
- $10^5$  to  $10^8$  viable cells or spores – cause illness
- $>10^3$  cells/g food – not safe for consumption



*Infections / outbreaks due to *B. cereus* are underreported because of usually mild disease (lasts < 24 h)*

# ISO 7932

## Colony count of *Bacillus cereus* at 30°C



Day  
1

25 g/ml TEST SAMPLE IN 225 ml BPW  
MAX. RECOVERY MEDIUM,  
RINGERS OR PEPTONE SALINE



0.1 ml  
(MYP Agar )  
30 °C FOR 18-48h



Day  
2,3

PICK OFF 5 CHARACTERISTIC COLONIES  
AND SUBCULTURE ON NUTRIENT AGAR  
35 OR 37 °C FOR 18-24 h

Day  
3,4

BIOCHEMICAL CONFIRMATION:  
GLUCOSE FERMENTATION (+), VOGES PROSKAUER (ACETOIN +),  
NITRATE (+)



# **B. cereus by Duopath® Cereus Enterotoxins** **Emetic B. cereus by Singlepath Emetic TOX MRK**

Day  
1

suspend  
1 - 3 presumptive colonies into 1  
ml CGY Broth + 1% Glucose,


4 h at 37°C

150 µl

**NO**

***B. cereus*  
not present**

**Emetic *B. cereus*  
not present**



**YES**

***B. cereus*  
present**

**Emetic *B. cereus*  
present**



**M.Y.P. Agar**



***Clostridium perfringens***

# INDICATOR ORGANISMS

## Clostridium perfringens

- In human faeces
- Sulphite reducing Clostridium
- Gram-positive, rod-shaped
- Anaerobic spore-forming bacteria
- Spores that are resistant to environmental conditions (temperature, pH, UV, water treatment processes, disinfection)
- Occurs smaller numbers than *E. coli*.



# ***C. perfringens***



- ***C. perfringens* is an indicator for fecal contamination**
- ***C. perfringens* spores indicate effectiveness of water purification processes**
- **Possible index parameter for occurrence of persistent pathogens e.g. viruses and oocysts (cryptosporidia)**

# ISO 7937:

## Colony count of *Clostridium perfringens*

Dilute food sample and add  
1 ml of sample to 15 ml  
TSC agar

TSC-Agar, anaerob  
(37 °C, (20 ± 2) h)

Select plates with less than 150  
colonies with black or brown  
colonies

Confirmation test  
Thioglycollate medium  
LS medium  
Nitrate motility medium  
Lactose gelatine medium



# ISO 7937: Alternative Identification procedure

## Colony count of *Clostridium perfringens*





Dilute food sample and add  
1 ml of sample to 15 ml  
TSC agar

TSC-Agar, anaerob  
(37 °C, (20 ± 2) h)

Select plates with less than 150  
comonies with black or brown  
colonies

Confirmation test  
Fluorescence  
reading

TSC-agar base



TSC agar MERCK  
supplemented  
with *Clostridium*  
*Perfringens* Selective  
Supplement:  
Cycloserine & MUP



Coliforms and E. coli


- **The coliform group consists of several genera of bacteria belonging to the family of Enterobacteriaceae. The historical definition of this group has been based on the method used for detection of lactose fermentation .**



# Indicator and index concept



- Mossel (1982) defined the term 'marker organism' which refers to two different functions, index and indicator.
- 'Index organisms' are related, directly or indirectly, either to the health hazards or to the presence of pathogens.
- On the other hand, 'indicator organisms' are related only to the effects of treatment processes or control of water quality.



## Total Coliforms (Marker Organisms)



### Indicator Organisms

Coliforms

indicate insufficient hygiene  
during a process/production

### Index Organisms

*E.coli*

indicates a potential hazard  
to health due to fecal contamination

# Gas and acid production from lactose



- *Escherichia coli*

- *Klebsiella*

- *Enterobacter*

- *Citrobacter*

- Gas + Acid from Lactose at 37°C:  
**95 % of Coliforms**


- Gas + Acid from Lactose at 44°C:  
**ONLY 90 % of E.coli**

- Tryptophanase:  
**99 % of E.coli are Indole positive**

- **The presence of fecal coliforms in drinking water indicates that an urgent public health problem probably exists, since human pathogens of co-exist with fecal coliforms.**
- **Therefore each total coliform-positive sample has to be analyzed if it contains fecal coliforms.**

# Reference Method: ISO 4831

## General guidance for coliforms – MPN technique



# Reference Method: ISO 7251 & ISO 11866-1

## Horizontal method for E.coli - MPN



3 tubes of LST broth  
Single strenght

3 tubes of LST broth  
Double strength

Incubation:  
24 to 48 hours at 37°C +/- 1°C

Check tubes for gas formation

**MPN**  
**24-48 h**

Subculture in Peptone water  
24-48 h at 44 ± 1°C

Subculture in EC Broth  
24-48 h at 44 ± 1°C

Indole


Gas

Coliform

-

E.coli

+




-

+

**Confirmation**  
**24-48 h**

# Reference Method: ISO 11866-1


## E.coli – MPN-MUG Foods



# Chromocult<sup>®</sup> TBX Agar

**Chromogenic Culture Medium for the detection and enumeration of  $\beta$ -glucuronidase-positive *Escherichia coli* in foods.**

**The culture medium complies with the recommendations in ISO 16649**



**Chromocult<sup>®</sup> TBX Agar # 1.16122**



# Chromocult® TBX Agar

## Chromocult® TBX Agar -

### New EU Food Regulation Law 2073/2074/2075/2076

**Chromocult® TBX Agar is used in the regulation for the detection of :**

- Living shells, snails
- Cooked crabs
- Fruits and vegetables
- Minced meat, Meat preparations
- Milk from cheese thermally treated
- Butter and cream not pasteurized



# Chromocult® TBX Agar

## Mineral - Modified Glutamate Agar (MMGA)

Culture medium for the resuscitation of severely stressed cells (for example by heating, drying, chemical preservation, freeze-drying, acidification) of *E. coli* from food samples (ISO 16649 - Part 1)

In combination with  
**Chromocult® TBX Agar**,  
the enumeration of damaged  
*E. coli* is possible.



**MMGA Agar # 1.09045**

**MMGA overlayed with a  
cellulose membrane**


# Chromocult<sup>®</sup> TBX Agar

## *E. coli* Enumeration in Processed Foods (ISO 16649-1) Colony-count Technique using Membranes

**Resuscitation**

**Selective medium**

**Counting**




# Chromocult<sup>®</sup> TBX Agar

## *E. coli* Enumeration in Fresh Foods (ISO 16649-2) Colony-count Technique using Pour Plate Method

**Selective medium**

**Counting**



# Reference method : ISO 4832

## General guidance for coliforms – colony count



**2 plates of VRB agar**  
**( 1 ml sample + 15 ml agar )**  
24 h at 30° or 35-37 °C

**Counting**  
**Confirmation**  
**24 h**



**Coliforms**


**red colonies , diameter > 0,5 mm**  
**reddish zone of precipitated bile**

# Rapid Methods for Water and Food


## Approved Alternative Methods Chromocult Coliform Agar

**USEPA  
for water**

**ISO 9308  
for water**




**AOAC-RI  
for food**



# Chromocult® Coliform Agar



Chromogenic Agar for the simultaneous detection of total coliforms and *E. coli* in drinking water and processed food samples.



**Chromocult® Coliform Agar # 1.10426**

# Chromocult® Coliform Agar


**Water Testing: Approved Method (USEPA/ ISO 9308**






# Chromocult® Coliform Agar

**Food Testing: AOAC™ Validated Method (Processed Food)**



# Chromocult® Coliform Agar ES

Chromogenic Agar for the simultaneous detection and colony count of total coliforms and *E. coli* in food samples.



**Chromocult® Coliform Agar ES # 1.00850**

# Chromocult® Coliform Agar ES


## Food Testing

**2007: AOAC™ validated method for the detection of total coliforms and *E. coli* in fresh food (raw ground beef, raw ground chicken, and raw milk)**




# Chromocult® Coliform Agar ES


## Food Testing: AOAC™ Validated Method (Fresh Food)



## Simplified Microbial Water Testing

A large, stylized black scroll that frames the central text. It starts at the top left, curves around the top, and ends at the bottom right where it meets a red ribbon.

# Readycult COLIFORMS



# Chromogenic Substrate Coliform Test

- **The chromogenic substrate coliform test is recommended for the analysis of drinking and fresh source water samples acc. to U.S. Standard Methods , 19th Edition 1995**
- **Readycult Coliforms is US-EPA approved since October 2002**



## Readycult® Coliforms Procedure

- Add one Readycult® snap pack to 100 ml of water in a sterile 120 ml to 150 ml bottle.
- Seal vessel and shake completely to dissolve granules.
- Incubate 24 hr +/- 1 hr at 35-37°C
- Interpretation of results



# Interpretation of Results

## NEGATIVE RESULTS (no coliforms)

**Broth Stays Clear After Incubation for 24 Hours**



**NEGATIVE**

**If Broth Becomes Turbid After Incubation the Sample is Still Negative. (Possible Growth of Non-coliforms)**



**NEGATIVE**



# Interpretation of Results

## POSITIVE TOTALCOLIFORMS

- Coliforms Present
- $\beta$ -D-Galactosidase Positive

**Distinct Color Change  
to Blue-Green**




# Interpretation of Results

## ***E. coli* POSITIVE**

If the Broth is Total Coliform Positive,  
(Blue-Green in Color) Check with UV  
Light for FLUORESCENCE.

If the Broth Fluoresces then:

- ***E. coli* is Present**
- $\beta$ -D-Glucuronidase Positive



# Interpretation of Results

*E. coli*

## BIOCHEMICAL CONFIRMATION WITH THE INDOLE REACTION

- Indole Reaction 99% Accurate.
- Tryptophanase Positive.

Overlay the Total Coliform or *E. coli* Positive Sample with Kovac's Indole Reagent. The Formation of a Red Ring Indicates a Positive Indole Test

INDOLE POSITIVE

**E. Coli Confirmed 99%**



# Enterococci

## Index Organisms



Intestinal Enterococci like *E. faecalis*, *E. faecium*, *E. durans* and *E. hirae* indicate a potential hazard to health due to fecal contamination




Lower concentration than *E. coli* in faeces



Survive longer than *E. coli* in environment

# Reference Method: Slanetz and Bartley Agar (ISO 7899-2)



**Counting  
Confirmation  
48 h**

## SLANETZ and BARTLEY Agar # 1.05262



*Enterococcus faecalis*

**Red** colonies

# Bile Aesculin Azide Agar # 1.00072

*Enterococcus faecalis*




**Red** colonies with black stained medium



**Black** color in the medium below the membrane

# Chromocult® Enterococci Agar

Chromogenic Culture Medium for the Detection and Enumeration of Enterococci in Water.



*Enterococcus faecalis*  
**RED** colonies

**Chromocult® Enterococci Agar # 1.00950**



# Chromocult® Enterococci Agar


## Water Testing

**2006: Accepted and approved  
as alternative method according to  
the EU Directive on Drinking Water in Germany  
Validated method ISO 17994**



# Chromocult® Enterococci Agar

## Water Testing: Approved Method (EU Directive)





# Rapid Detection of *E. coli* O157 and EHEC's



# EHEC's: Emerging Pathogens

- Variety of *E. coli* present in nature, usually found in the intestines of healthy humans and animals and offer beneficial properties.
- Highly pathogenic *E. coli* (EHEC) strains exist.
- *E. coli* O157:H7 is most prominent EHEC strain which can cause severe foodborne illness in a VERY low dosage (10 CFU can cause bloody diarrhea and abdominal cramps).
- Complication: The elderly and children <5 years are highly susceptible to Hemolytic Uremic Syndrome (HUS: destruction of red blood cells and kidneys failure).
- Disease is mainly caused by production of Verotoxins (= Shiga-like toxins)



# EHEC - Relevance to Public


ca. 100,000 people in the US are annually affected

- Fecal contaminated meats (most frequent outbreaks)
- undercooked beef (HAMBURGERS!!!!!!)
- water
- raw milk
- unpasteurized apple juice/cider
- sandwiches
- lettuce
- dry cured salami
- produce from manure-fertilized gardens
- radish sprouts, alfalfa sprouts




**EHEC disease**  
**also known**  
**as *Hamburger***  
***Disease***



# ISO Standard for testing of *E. coli* O157 in Food



# Singlepath® *E. coli* O157 for the Detection of *E. coli* O157 in Food

Day 1	25 g/ml TEST SAMPLE 225 ml in <b>m-TSB</b> 37 °C FOR 8 h		OR	25 g/ml TEST SAMPLE in 225 ml <b>mEC+N or mTSB+N</b> 37- 42 °C FOR 24 h	
Day 2	<b>NEW</b> <b>No</b> <i>E. coli</i> O157 present  <b>only 8 h</b>				<b>YES</b> <i>E. coli</i> O157 present  <b>only 1 medium</b>
Day 3	If positive, confirmation by culture and biochemical tests required			<b>STREAK OUT ON</b> <b>CT-SMAC Agar</b> 37 °C FOR 24 h	
					

# Detection of E. coli O157:H7 in drinking water

**NEW**



The first 24 hours Rapid Method for the detection of E. coli O157:H7 in water using ReadyCult Coliforms for screening and Singlepath E. coli O157 for confirmation



# Rapid Detection of *E. coli* O157 in WATER using Singlepath® *E. coli* O157

**100 ml water sample  
+ Readycult Coliforms 100  
Enrichment 18-24 h, 35-37°C**



**Blue-Green colour  
Fluorescence negative  
*E. coli* O157:H7/:H- may  
be present**



**Red ring formed (Indole  
pos.)  
*E. coli* O157:H7/:H-  
probably present**



**Blue-green colour  
No Fluorescence  
Indol positive**




**150 µl of Indole-pos. sample  
immediately transferred to Singlepath  
*E. coli* O157;  
<20 min RT  
*E. coli* O157:H7/:H- is present**

**YES  
*E. coli* O157  
present !**

# Duopath® Verotoxins

## Unique Lateral Flow Test for Verotoxin Identification


- Confirmation test for the detection of Verotoxins (Shiga-like Toxins) 1 and 2 from pathogenic *E. coli* isolates (VTEC's or STEC's)
- 6 h enrichment in CAYE Broth required
- No centrifugation step
- No boiling
- 10 min. treatment with Polymixin recommended to free cell-associated Verotoxins
- read results within 20 min  
(most results within 5-10)
- Lower detection limit for toxins:  $10^7$  cfu/ml
- Approvals: AOAC-RI, CE, FDA 510(k), IVD (Japan)




# Confirmation of VTEC's using Duopath® Verotoxins

day  
1

**SMAC Agar**



**CT-SMAC Agar**



Streak out on BHI Agar + Carb.  
at 37°C for 24 h, suspend in  
250µl PBS, Polymyxin 10 min

**OR**

Suspend 1 – 5 colonies  
into CAYE Broth +  
Carbadox at 37°C for 6 h,  
Polymyxin 10 min



**NO  
Verotoxin  
present**



**YES  
Verotoxins detected  
and presence of  
VTEC  
confirmed !**

# Duopath® Legionella



## Natural habitat of Legionella

Legionella are frequently found in aquatic environments and some species have been recovered from soil.


Temperature is a critical determinant for Legionella proliferation.

Colonization of hot water tanks is more likely if tank temperatures are between 40 and 50°C.



# Detection and enumeration of Legionella

## Direct membrane filtration method (ISO 11731-2)



# Confirmation of presumptive Legionella colonies (ISO 11731-2)



BCYE / GVPC Agar

**Presumptive Legionella**  
Typical colonies on BCYE / GVPC Agar



**BCYE with Cystein**



**BCYE without Cystein**



**Growth**

48 h / 36 ± 2°C



**No growth**



**Unsafe Method as some Legionellae can grow on BCYE without Cys:  
Leads to false-negative results**



# Thank you for your attention

